

# CATH•SECURE® Line of Products



## Features and Benefits

- Secures different types of tubing and sizes
- Improves patient's quality of life
- Quick and Easy to remove
- Cost Effective • SKU Reduction
- Increases nursing productivity
- Hypoallergenic

## Product Info

- **Item No. 5445-2**  
The Original Multi-Purpose and Foley Catheter Securement Device
- **Item No. 5445-3**  
Multi-Purpose and Foley Catheter Securement Device with Extended Tab
- **Item No. 5445-4**  
Wraps Around Technology to Limit Pistoning
- **Item No. 5445-6**  
Water-Resistant Multi-Purpose Tube Holder

**DEHP FREE • LATEX FREE**

Distributed by Pat Coakley Medical  
Email: [sales@patcoakleymedical.ie](mailto:sales@patcoakleymedical.ie)



Units E3 & E4  
Network Enterprise Park  
Kilcoole, Co. Wicklow

Tel: +353 1 2018070  
Fax: +353 1 2018074

5445-2



CATH-SECURE is the original multi-purpose securement device, providing an efficient, standardized system for securing catheters and medical tubing.

Package: 50 Box / 8 Case  
2.5" Long Tab & 3" x 1.75" Wide Base  
DEHP Free • Hypoallergenic  
Latex Free

GMS Code: 93910

5445-3



The CATH-SECURE multi-purpose tube holder features an extended tab for larger size tubing and drainage ports while providing a comfortable, cost effective solution for your securement needs.

Package: 50 Box / 8 Case  
2.75" Long Tab & 3" x 1.25" Wide Base  
DEHP Free • Hypoallergenic  
Latex Free

GMS Code: 93909

5445-4




The CATH-SECURE DUAL TAB extends from the center and wraps around the tubing to prevent pistoning. The CATH-SECURE DUAL TAB is lightweight and flexible, improving the patient's comfort and mobility.

Package: 50 Box / 8 Case  
(2 Dual) 3" Wrap Tabs x 3" Wide Base  
DEHP Free • Hypoallergenic  
Latex Free

GMS Code: 93907

5445-6



The CATH-SECURE PLUS, our newest patented securement device, is made using a breathable, water-resistant, hypoallergenic material that features a butterfly design, improving overall clinical outcomes.

Package: 50 Box / 8 Case  
2.5" Long Tab & Butterfly Base  
DEHP Free • Hypoallergenic  
Latex Free

GMS Code: 93915

7334-1



**TUBE-SECURE** *for Kids*<sup>TM</sup>

The TUBE SECURE for Kids provides an efficient, standardized system for securing pediatric medical tubing.

Package: 50 Box / 8 Case  
1.87" Base, 1.31" x .25" Tab  
DEHP Free • Hypoallergenic  
Latex Free

GMS Code: 93913

7334-2



**CATH-SECURE**<sup>®</sup> *for Kids*

The CATH-SECURE<sup>®</sup> for Kids is a multi-purpose tube holder that provides an efficient, standardized system for securing pediatric catheters and medical tubing.

Package: 50 Box / 8 Case  
1.75" Long Tab & 1.5" x .75" Wide Base  
DEHP Free • Hypoallergenic  
Latex Free

GMS Code: 94031

