



Guedel Airways

Description

A single-use oropharyngeal airway, manufactured from LLDPE (LINEAR LOW-DENSITY POLYETHYLENE) to the highest British Specifications conforming to: B.S. EN 12181 1998. Design incorporates a colour coded bite block for ease of selection of colour coded catheters.



Sizes

Code	Sizes
38/000	000 None
38/00	00 Blue
38/0	0 Black
38/1	1 White
38/2	2 Green
38/3	3 Orange
38/4	4 Red

Indication

Designed to maintain a patent oropharyngeal airway for patients under general anaesthetic during administration of the anaesthetic, or during recovery.

Advantages

- The smooth finish of the airway reduces the risk of damage to oropharyngeal tissues.
- Bite block reinforcement prevents occlusion of airway in the event of involuntary biting.
- Available in seven sizes.
- All our product range is latex and phthalate free.

Presentation

- Supplied sterile in individual single wrapped pouches, 50 pouches per carton.

Sterilization

- Sterilized by gamma irradiation, single use only.