


NEW!


EASY CLEANING BASE

STAINLESS STEEL HOOKS



STAINLESS STEEL
LOCKING SLEEVE

OPTIONAL FUNNY CASTORS FOR CHILDREN



EASY CLEANING



Description

Serum infusion stand with plastic base on wheels (500mm width base) easy to clean, without harshness or cavity. Available in 6 colors, black, grey, blue, red, yellow and green. Lower tube and rod telescopic sleeve $\varnothing 25\text{mm}$ $\varnothing 16\text{mm}$ 18/10 stainless steel AISI 304. Locking of the telescopic rod with brass screw on robust steel sleeve. Maximal capacity load 2 kg per hook.

| Ref. | Description | Height maxi |
|---------|------------------------------------|-------------|
| 1773 P | With 2 Polyamid hooks | 2000mm |
| 1773 I | With 2 stainless steel hooks | 2000mm |
| 1773 4 | With 4 stainless steel hooks | 2000mm |
| 1773 U | With 2 stainless steel hooks « U » | 2000mm |
| 1773 4U | With 4 stainless steel hooks « U » | 2000mm |

Choose your base color by using a letter as following :

B for **Blue** / G for **Grey** / R for **Red** / J pour **Yellow** / V pour **Green**




RANGE OF INFUSION STAND

■ 5 TYPES OF STANDS AVAILABLE

Cast iron stand
Covered with 18/10 stainless steel
Ø 280 mm



4 castors stainless steel stand
Square tube 35 x 35 mm
600 mm width base



5 castors stainless steel stand
Square tube 25 x 25 mm
600 mm width base



5 castors moulded Nylon stand
520 mm width base



5 castors aluminium stand
Polished moulded aluminium
520 mm width base



■ CASTORS FOR STAINLESS STEEL STANDS

2 types of castors are available



Twin black castors
Ø 50 mm (authorized load
40 kg/castor)



Antistatic castors
Ø 50 mm Subject to regulation with regard to
gas anaesthetics in ops room

- Weighting available (screwed on the bottom of stand). Obligatory for infusion stands with a maximum height of 2400 mm.

■ CASTORS FOR NYLON AND ALUMINIUM STANDS

5 types of castors are available



Twin black castors Ø 50 mm
(authorized load 40 kg/castor)



Twin black castors Ø 50 mm (authorized load 40 kg/castor)
Associated with coloured hooks and sleeves.

- Weighting available (screwed on the bottom of stand). Obligatory for infusion stands with a maximum height of 2400 mm.

■ STANDARD SLEEVE

For all stands


18/10 stainless steel tube Ø 30 mm –
Length 1000 mm
With stainless steel locking sleeve.



■ SERUM RODS

18/10 stainless steel Ø 16 mm tube
1000 mm length allowing height adjustment to 2000 mm
1500 mm length allowing height adjustment to 2400 mm

5 types of safety hooks



2 unbreakable polyamide hook
3 kg authorized per hook



2 stainless steel hooks
authorized load 12 kg
(6 per hook)



4 stainless steel hooks
authorized load 12 kg
(3 per hook)



4 plastic hooks
authorized load 12 kg
(3 per hook)

■ ACCESSORIES

To meet the requirements of the different departments a wide range of accessories is available at the end of this brochure (instrument trays, serum holders...)



Réf. 770 I



Réf. 770 4

2 STAINLESS STEEL SAFETY HOOKS

| Reference | Maxi height |
|---------------|-------------|
| 770 I | 2000 mm |
| 770 IL | 2400 mm |

4 STAINLESS STEEL SAFETY HOOKS

| Reference | Maxi height |
|---------------|-------------|
| 770 4 | 2000 mm |
| 770 4L | 2400 mm |



Réf. 1774 I



Réf. 1774 P

2 POLYAMIDE SAFETY HOOKS

| Reference | Maxi height |
|---------------------------------------|-------------|
| 1774 P with twin castors | 2000 mm |
| 2774 P with antistatic castors | 2000 mm |

2 STAINLESS STEEL SAFETY HOOKS

| Reference | Maxi height |
|---------------------------------------|-------------|
| 1774 I with twin castors | 2000 mm |
| 2774 I with antistatic castors | 2000 mm |

4 STAINLESS STEEL SAFETY HOOKS

| Reference : | Maxi height |
|---------------------------------------|-------------|
| 1774 4 with twin castors | 2000 mm |
| 2774 4 with antistatic castors | 2000 mm |



Twin black castors Ø 50 mm or
Antistatic castors Ø 50 mm
Black color



Ref. 2775 I

Ref. 2775 4

Twin black castors Ø 50 mm or
Antistatic castors Ø 50 mm
Black color



Ref. 1775 P



■ 2 POLYAMIDE SAFETY HOOKS

| Reference | Maxi height | Weighted |
|-----------|-------------|----------|
|-----------|-------------|----------|

With twin castors

| | | |
|----------------|---------|-----|
| 1775 P | 2000 mm | no |
| 3775 P | 2000 mm | yes |
| 3775 PL | 2400 mm | yes |

With antistatic castors

| | | |
|---------------|---------|-----|
| 4775 P | 2000 mm | yes |
|---------------|---------|-----|

■ 2 STAINLESS STEEL SAFETY HOOKS

| Reference | Maxi height | Weighted |
|-----------|-------------|----------|
|-----------|-------------|----------|

With twin castors

| | | |
|----------------|---------|-----|
| 1775 I | 2000 mm | no |
| 3775 I | 2000 mm | yes |
| 3775 IL | 2400 mm | yes |

With antistatic castors

| | | |
|---------------|---------|-----|
| 2775 I | 2000 mm | no |
| 4775 I | 2000 mm | yes |

■ 4 STAINLESS STEEL SAFETY HOOKS

| Reference | Maxi height | Weighted |
|-----------|-------------|----------|
|-----------|-------------|----------|

With twin castors

| | | |
|----------------|---------|-----|
| 1775 4 | 2000 mm | no |
| 3775 4 | 2000 mm | yes |
| 3775 4L | 2400 mm | yes |

With antistatic castors

| | | |
|----------------|---------|-----|
| 2775 4 | 2000 mm | no |
| 4775 4 | 2000 mm | yes |
| 4775 4L | 2400 mm | yes |

Distributed by

Pat Coakley Medical Ltd

sales@patcoakleymedical.ie

www.patcoakleymedical.ie

phone: 01 2018070 / 2



STAINLESS STEEL STAND 5 CASTORS