


Great Bear[®] Healthcare

A UK manufacturer committed to Urology



Great Bear Product Guide

Hello...
We are Great Bear.

Contents

Libra Leg Bags

Libra Conform Leg Bags

Paediatric Leg Bags

Single Use Night Bags

Reusable Night Bags

LibraSheath

GB Soft Skin[®]

Libra Lever Catheter Valve

Libra Leg Bag Straps

Importance Of Using A Fixation Device

Libra Leg Bag Holders

Fix-it Retaining Strap

CubCath Intermittent Catheter

Libra Leg Bag

Utilising high quality materials with a user friendly design, the Libra Leg Bag provides total peace of mind for users and carers alike.

With a comprehensive range of sizes and tube lengths, we are confident that there will be a Libra Leg Bag to suit all needs.

Libra Leg Bag – Prescription

All Prescription Libra Leg Bags come with Fabric Backing and one pair of Leg Bag Straps as standard


Product Code	Capacity	Tube Length	Quantity	PIP Code
10130G	350ml	Direct Inlet	10	
10131H	500ml	Direct Inlet	10	
10132J	750ml	Direct Inlet	10	
10133K	1000ml	Direct Inlet	10	
10100A	350ml	Short (10cm)	10	
10102C	500ml	Short (10cm)	10	
10106G	750ml	Short (10cm)	10	
10104E	500ml	Long (30cm)	10	
10108J	750ml	Long (30cm)	10	
10120W	1000ml	Long (30cm)	10	
10105F	500ml	Adjustable (45cm)	10	
10110L	750ml	Adjustable (45cm)	10	
10122Y	1000ml	Adjustable (45cm)	10	

Libra Leg Bag – Hospital

All Hospital Libra Leg Bags come **without** Fabric Backing and one pair of Leg Bag Straps in every peel pouch as standard.

Product Code	Capacity	Tube Length	Quantity
10230C	350ml	Direct Inlet	10
10200W	350ml	Short (10cm)	10
10202Y	500ml	Short (10cm)	10
10206C	750ml	Short (10cm)	10
10204A	500ml	Long (30cm)	10
10208E	750ml	Long (30cm)	10
10205B	500ml	Adjustable (45cm)	10
10210G	750ml	Adjustable (45cm)	10

Key Features



Did you know... The Libra Leg Bag tap can be swivelled for comfort and extra security. It also features raised ridges to allow the visually impaired to identify whether the bag is open or closed.

Libra Conform Leg Bag

Designed for maximum discretion and comfort, the Libra Conform is a three chambered Leg Bag that will wrap around the contour of the leg.

Libra Conform – Prescription

All Prescription Libra Conform come with Fabric Backing and one pair of Leg Bag Straps as standard

Product Code	Capacity	Tube Length	Quantity	PIP Code
10114Q	500ml	Short (10cm)	10	
10115R	500ml	Long (30cm)	10	
10116S	750ml	Short (10cm)	10	
10117T	750ml	Long (30cm)	10	

“The Libra Conform fits perfectly under my jeans – finally a product that works”

Patient, Nottinghamshire

Key Features

- The three chambers enhance comfort and discretion as the “ballooning” effect is reduced – less noticeable under clothing.
- Urine is held across the three chambers, minimising any audible noises of the urine inside the bag.
- Weight of urine is more evenly distributed – making the Libra Conform more comfortable to wear.



Did you know... The Libra Conform shares all the great features of the Libra Leg Bag – including its needle-free sampling port, fabric backing and easy to use lever action tap.

Paediatric Leg Bag

Having to use a Urinary Drainage Bag can be an overwhelming and daunting experience – least of all when you are only young. Great Bear's Paediatric range is designed specifically for babies and young children with a fun child-friendly appearance to make the experience less traumatic.

Paediatric Leg Bags- Prescription

All Prescription Paediatric Bags come with Fabric Backing and one Leg Bag Strap as standard

Product Code	Capacity	Tube Length	Quantity	PIP Code
10101B	150ml	Direct Inlet	10	
10125B	150ml	Adjustable (45cm)	10	
10111M	250ml	Direct Inlet	10	
10135M	250ml	Adjustable (45cm)	10	
10103D	350ml	Direct Inlet	10	
10151D	350ml	Adjustable (45cm)	10	



Please Note!


Great Bear's Paediatric range comes with their own specific Leg Bag Strap. One strap is included in the box, but they are also available to order in packs of five.



Paediatric Straps- Prescription

Product Code	Size	Quantity	PIP Code
10633P	30cm	5	
10657P	50cm	5	

Key Features



Single Use Night Bags

Our Night Drainage Bags can be connected to a Leg Bag to allow overnight drainage. The single use design means that, once emptied, they can be disposed of for ease and convenience.

Single Use Night Bags – Prescription

Code	Name	Capacity	Tube Length	Description	Qty.	PIP Code
10300S	GB2	2-Litre	90cm	Non Sterile and Non Drainable	10	
10303V	GB3	2-Litre	90cm	Non Sterile with Single Use T-Tap	10	
10GB3S	GB3S	2-Litre	90cm	Sterile with Single Use T-Tap	10	
10330Y	GB5	3-Litre	120cm	Non Sterile with Single Use T-Tap	10	

Single Use Night Bags – Hospital

Code	Name	Capacity	Tube Length	Description	Qty.
10301T	GB2	2-Litre	90cm	Non Sterile and Non Drainable	10
10304W	GB3	2-Litre	90cm	Non Sterile with Single Use T-Tap	10
10GB3H	GB3S	2-Litre	90cm	Sterile with Single Use T-Tap	10



Please Note!

All single use Night Bags must be connected to a Leg Bag for overnight drainage. **DONOT** connect directly to a Catheter



Key Features

- Drainable and non drainable options available
- Choice of 2-litre and 3-litre capacities
- Can choose between sterile and non-sterile designs
- Punched out hanger holes to allow positioning on a Night Stand (see page 10)
- Tear open strip on the GB2 and Single Use T-Tap on the GB3, GB3S and GB5 to allow for quick, easy emptying
- **TamperProof** – Our single use Night Bag range **cannot** be reused after emptying
- A safety clip on the GB3, GB3S and GB5 prevent the tap from being accidentally activated prior to emptying



Fact... 80% of all UTIs in hospitals are attributed to Indwelling Catheters¹. Single Use Night Bags might be the safest way to reduce the risk of infection²

¹NHS Institute for Innovation & Improvement; ²Jones S, Brooks A, Foxley S, Dunkin J (2008) "Evidence: Care of urinary catheters & drainage systems" Nursing Times 103:42,48-50.



Did you know... Our GB3S is the GB3 Night Bag in sterilised form. It is designed specifically for those who are more vulnerable and susceptible to infection but who like the convenience of a Single Use Night Bag

Reusable Night Bags

Unlike the single use options, the GB4, GB4 Fold Up and GB6 can be reused over a period of 5–7 days. They can either be connected to a Leg Bag as part of an overnight linksystem or they can be connected directly to the Catheter or Sheath.

Reusable Night Bags – Prescription

Code	Name	Capacity	Tube Length	Description	Qty.	PIP Code
10400W	GB4	2–Litre	100cm	Sterile and reusable with Lever Action Tap	10	
10402Y	GB4 Fold Up 2–Litre	2–Litre	100cm	Sterile and reusable with star tubing and foldable Lever Action Tap	10	
10430U	GB6	3–Litre	120cm	Sterile and Reusable with star tubing and Lever Action Tap	10	

Reusable Night Bags – Hospital

Code	Name	Capacity	Tube Length	Description	Qty.
10401X	GB4	2–Litre	100cm	Sterile and reusable with Lever Action Tap	10
10403Z	GB4 Fold Up	2–Litre	100cm	Sterile and reusable with star tubing and foldable Lever Action Tap	10

Night Stands

Code	Product Name
10708S	3 in 1 Metal Night Stand
10705P	Great Bear Flat Pack Hanger
10700J	Great Bear Plastic Bed Hanger (Hospital Only)

Please Note!

It is important to use a Night Stand for effective drainage overnight. Great Bear have both metal and flat pack hangers available to order. Night Stands are provided to patients complimentary upon registration with our Home Delivery Service.



GB4 Fold Up

GB4



In closer detail... The tubing on the GB4 Fold Up and GB6 has an integral star formed design. This ensures free flow of urine, even when the tubing is kinked. The outer surface remains smooth – harbouring no nasty bacteria and being gentle on the skin.



Did you know... The tap on the GB4 Fold Up can be folded up to prevent it having prolonged contact with surfaces that may be contaminated.

Libra Sheath

The Libra Sheath is available in a range of styles and sizes and is designed to help men manage their continence issues.

Libra Sheath – Prescription

All Libra Sheath come with a measuring guide and 30 hair protectors as standard.

Code	Style	Size	Quantity	PIP Code
GBLSS24	Standard	24mm	30	
GBLSS28	Standard	28mm	30	
GBLSS31	Standard	31mm	30	
GBLSS33	Standard	33mm	30	
GBLSS35	Standard	35mm	30	
GBLSS40	Standard	40mm	30	
GBLSW24	Wide Band	24mm	30	
GBLSW28	Wide Band	28mm	30	
GBLSW31	Wide Band	31mm	30	
GBLSW33	Wide Band	33mm	30	
GBLSW35	Wide Band	35mm	30	
GBLSW40	Wide Band	40mm	30	
GBLSP24	Pop On	24mm	30	
GBLSP28	Pop On	28mm	30	
GBLSP31	Pop On	31mm	30	
GBLSP33	Pop On	33mm	30	
GBLSP35	Pop On	35mm	30	
GBLSP40	Pop On	40mm	30	

Selecting your Libra Sheath

Style 1: Standard



Style 1: Standard

A standard length Sheath with centrally positioned adhesive. Libra Sheath Standard are suitable for men with an average length penis.

Style 2: Pop On



Style 2: Pop On

A shorter length Sheath combined with a maximum adhesive area will provide men with a short or retracted penis a suitable Sheath option.

Style 3: Wide




Style 3: Wide Band

A standard length Sheath with a greater adhesive area to maximise security – particularly for those who are more active and have trouble keeping other Sheaths in place.

Key Features

- Integrated adhesive eliminates the need for adhesive strips and tape
- The nozzle is kink resistant to maintain continual flow of urine into the Drainage Bag
- Silicone is malleable to provide a secure and comfortable fit to the penis
- Manufactured using high quality silicone that is breathable – maximum comfort and skin protection
- 30 hair protectors and a measuring guide in each box as standard
- 100% Latex Free



“They’re amazing!”

Company Nurse, Forest Of Dean

GB Soft Skin®

The GB Soft Skin® range consists of silicone based Medical Adhesive Removers and Barrier Wipes that are designed to be used in conjunction with a Medical Adhesive product – such as a Sheath or Stoma Pouch.

The Barrier Wipes are to be used prior to application of a Medical Adhesive product, whilst an Adhesive Remover allows for that product to be gently removed afterwards – leaving no sticky residue behind.

GB Soft Skin Range® – Prescription

Code	Description	Quantity	PIP Code
10380A	GB Soft Skin® Medical Adhesive Remover Spray	50ml Aerosol	
10381B	GB Soft Skin® Barrier Wipes	30	
10382C	GB Soft Skin® Medical Adhesive Remover Wipes	30	



Did you know... The Adhesive Remover wipe is really good for removing any built up residue on the skin – ideal for Stoma Pouches or even wound dressings.

Key Features – Barrier Wipes

- Alcohol-Free, non-sting formula
- Barrier Wipes leave a protective layer on the skin. This provides a “barrier” between the skin and Medical Adhesive product to prevent irritation
- Maximises the effectiveness of the product in use
- Non-drip – dries quickly to allow the Sheath or Stoma Pouch to be applied instantly



Key Features – Adhesive Remover

- Available in either an aerosol spray or individual wipe sachets
- Alcohol-Free, non-sting formula
- Removal of adhesives can cause trauma to the skin – GB Soft Skin® alleviates this discomfort to ensure a pain-free removal
- The spray is non-chill, dries quickly and leaves no white powdery residue
- The aerosol can be sprayed upside down to allow access to hard-to-reach places

Libra Lever Catheter Valve

With its easy to use opening and closing system, the Libra Lever Catheter Valve provides patients a secure, more comfortable and discreet alternative to Drainage Bag, where there is adequate bladder tone, capacity and sensation².

This sterile, tap-like device has been designed to fit all urethral and suprapubic Catheters and Sheaths.


²National Institute for Health and Clinical Excellence (NICE) 2012.

Libra Lever Catheter Valve – Prescription

Code	Description	Quantity	PIP Code
10540A	Libra Lever Catheter Valve – Lever Action Tap	5	

Libra Lever Catheter Valve – Hospital

Code	Description	Quantity
10541B	Libra Lever Catheter Valve – Lever Action Tap	5



Libra Leg Bag Straps

Libra Leg Bag Straps secure a Leg Bag to the leg. A Velcro style fastening allows the strap, and the Leg Bag, to be held in place.

Libra Leg Bag Straps – Prescription

Code	Size	Quantity	PIP Code
10622C	70cm:30cm – for use on the calf	5 pairs	
10637T	60cm:40cm – for use on the thigh	5 pairs	



Top Tip. Not sure on size? As a rule of thumb use a calf strap with a Long Tube or Adjustable Tube Leg Bag and a thigh strap with a Direct Inlet or Short Tube.



Key Features

- The two wavy silicone lines on the reverse of the strap will assist in securing the Leg Bag to the leg
- 25mm wide with a tapered Velcro fastening – designed for maximum comfort
- Tapered ends will not dig into the skin and allows for easy fastening and unfastening
- Stretchable – a comfortable and flexible fit
- Latex Free

Please Note!

Lack of Catheter securement can lead to damage and inflammation to the Urethra and Meatus. This could cause pain, discomfort and a high risk of infection for the patient³

³ Fisher, J. (2010)
"The importance of effective catheter securement"
British Journal of Nursing, 19:8



Libra Leg Bag Holder

Another popular method of securing a Leg Bag is the Libra Leg Bag Holder. With its 30% cotton count it will comfortably support Leg Bags whilst allowing the skin to breathe.

Libra Leg Bag Holder

Code	Colour Band	Size	Quantity	PIP Code
10660S	Yellow	Small (35–45cm)	4	
10655M	Blue	Medium (45–65cm)	4	
10654L	Brown	Large (55–80cm)	4	
10665X	Green	Extra Large (65–90cm)	4	



Top Tip... Measure the circumference of your thigh or calf at the widest point to find the correct size.



Key Features

- Encases the Leg Bag, offering maximum support, comfort and security
- Available in a range of sizes to ensure a good fit
- Colour coded for easy identification of size
- Has a higher cotton count than other Leg Bag Holders
– lessens skin irritation and allows the skin to breathe
- Universal design so can be used with any Leg Bag

Fix-it Retaining Strap

The Fix-it Retaining Strap is designed to secure Indwelling and Suprapubic Catheters, as well as Leg Bag tubing with a Sheath System.

Proper securement is essential to prevent urethral trauma.

Fix-it Retaining Strap

Code	Size	Quantity	PIP Code
10644A	Medium (45cm)	5	
10645B	Long (80cm)	5	
10646C	Short (35cm)	5	
10647D	Adjustable (150cm)	5	



"It is important to secure the Catheter after insertion to prevent movement & urethral traction"

European Association of Urology Nurses, 2012

Key Features

- Available in a range of sizes giving more choice
- It has a secure velcro fastening which is simple to use
- The 2 silicone lines of the reverse ensure a secure fit and eliminate slippage
- Raised tab on the silicone Fix-ittie allows for a better grip for the individual
- Prevents against pulling on the Catheter that could lead to pain and urethral trauma
- Clear labelling to denote size and washing instructions
- 100% Latex Free



Did you know... A user guide with easy to follow instructions and diagrams is available to ensure the Fix-it is used correctly. Email: marketing@gbhl.co.uk for a copy.


CubCath

A sterile, hydrophilic coated Catheter with smooth eyes and lubricious coating. The CubCath ensures comfort on insertion and removal.

CubCath – Prescription

Code	Length	Charriere Size	Funnel colour	Quantity	PIP Code
AVP06	Paediatric (28cm)	06ch	Green	30	
AVP08	Paediatric (28cm)	08ch	Blue	30	
AVP10	Paediatric (28cm)	10ch	Black	30	
AVF08	Female (18cm)	08ch	Blue	30	
AVF10	Female (18cm)	10ch	Black	30	
AVF12	Female (18cm)	12ch	White	30	
AVF14	Female (18cm)	14ch	Green	30	
AVM10	Male (40cm)	10ch	Black	30	
AVM12	Male (40cm)	12ch	White	30	
AVM14	Male (40cm)	14ch	Green	30	
AVM16	Male (40cm)	16ch	Orange	30	

Key Features



Great Bear Healthcare Ltd

1 Lambourne Crescent,
Cardiff Business Park,
Cardiff, CF14 5GF

www.patcoakleymedical.ie 01-2018070 sales@patcoakleymedical.ie

Sole Distributor in Ireland

Pat Coakley Medical Ltd

Units E3 E4
Network Enterprise Park
Kilcoole A63 H244
Co Wicklow
IRELAND