



PHCOM030 Commode Pan Liner

NHS Code: FSL297

- ✓ Pulp
- ✓ Recycled
- ✓ Disposable
- ✓ Single Use

Commode pan liner made of recycled and disposable pulp material with added waterproofing. High quality product designed for enhanced patient acceptance and comfort. Designed to fit all models of commode chairs and suitable for disposal in any macerator.

Benefits

- ✓ **Disposable product** for easy and quick handling, increased hygiene and reduced risk of contamination
- ✓ **Pulp product** made of recycled materials and disposable through macerators
- ✓ **Fits standard dispensing racks**
- ✓ **Controlled pouring** to ease handling
- ✓ **Cover lids** for added infection control and hygiene
- ✓ **Steady when standing** due to flat bottom surface
- ✓ **Strong sides** so able to support weight of contents when full
- ✓ **Easy to store** (product is stackable) and separate
- ✓ **Easy to carry** when full

Specification

- ✓ Capacity: 1700ml
- ✓ Material: Paper pulp (100% newsprint)
- ✓ Dimensions: 275L x 275W x 100H (mm)
- ✓ Weight: 42g ± 4g
- ✓ Colour: Grey
- ✓ Compliant with:
 - Medical Devices Directives 93/42/EEC Class 1 (2007/47/EC)
 - British standard PAS29:1999
- ✓ Disposability: By maceration

Fluid Testing


- ✓ Fluid Resistance: British Standard PAS29:1999
- ✓ Fluid Tested: Water
- ✓ Resistance Time >4 hours
- ✓ Resistance Temperature 35°C (±3°C)

Packaging



Dimensions: 550mm L x 275mm W x 785mm H

- Recyclable packaging
- Keep away from sunlight
- Manufacturer
- Product conforms with applicable EC directives
- Keep dry
- Do not re-use



Distributed by
 Pat Coakley Medical Ltd
 sales@patcoakleymedical.ie
 www.patcoakleymedical.ie
 phone: 01 2018070 / 2



Packing

Outer: 200 pieces per case
Pallet: 10 cases per pallet



BSI Kitemark for compliance with PAS29:1999

Barcodes



PHCOM030/APR2012



0044 (0) 208 8335 3636



0044 (0) 208 8335 0300



sales@hpchealthline.co.uk



HPC Healthline Limited, Colwood House, 158 Garth Road, Morden, Surrey, SM4 4LZ, England



www.hpchealthline.co.uk



www.facebook.com/hpchealthline



#hpchealthline

HPC Healthline