



PHTRA040

Medium Autoclavable Tray

NHS Code: FWK097

✓ Pulp ✓ Recycled ✓ Disposable ✓ Single Use ✓ Multipurpose

Autoclavable tray made of recycled and disposable pulp material with added waterproofing. High quality product designed for a multitude of functions such as carrying instruments, dressings and suitable for disposal in any macerator.

Benefits

- ✓ **Disposable product** for easy and quick handling, increased hygiene and reduced risk of contamination
- ✓ **Pulp product** made of recycled materials and disposable through macerators
- ✓ **'Autoclavable'** (by 3 sterilisation methods) without loss of form or function
- ✓ **Heat pressed** for a smoother finish and flatter surface, minimising residual particles
- ✓ **Steady when standing** due to flat bottom surface
- ✓ **Rolled edges** for improved strength when carried full
- ✓ **Easy to store** (product is stackable) and separate

Specification

- ✓ Material: Paper pulp (80% newsprint/20% virgin fibre)
- ✓ Dimensions: 181L x 137W x 20H (mm)
- ✓ Weight: 12.5g
- ✓ Colour: Grey
- ✓ Compliant with:
 - Medical Devices Directives 93/42/EEC Class 1 (2007/47/EC)
 - British standard PAS29:1999
- ✓ Disposability: By maceration

Fluid Testing

- ✓ Fluid Resistance: British Standard PAS29:1999
- ✓ Fluid Tested: Water
- ✓ Resistance Time >4 hours
- ✓ Resistance Temperature 35°C (±3°C)

Packaging



Dimensions: 800mm L x 385mm W x 415mm H



Recyclable packaging



Keep away from sunlight



Manufacturer



Product conforms with applicable EC directives



Keep dry



Do not re-use



Packing

Sleeve: 195 trays per sleeve
Outer: 1170 trays per case
Pallet: 9 cases per pallet
6 sleeves per carton

Barcodes



5 024951 010678 >

PHTRA040/JUN2012



0044 (0) 208 8335 3636



0044 (0) 208 8335 0300



sales@hpchealthline.co.uk



HPC Healthline Limited, Colwood House, 158 Garth Road, Morden, Surrey, SM4 4LZ, England



www.hpchealthline.co.uk



www.facebook.com/hpchealthline



#hpchealthline

HPC Healthline