







VITALCARE AUTO SAFETY LANCET

AUTO TYPE	NPC CODE	IMAGE DEPTH	PENETRATION	TYPE & SIZE	VOLUME OF BLOOD	PACKING	SUGGESTED TEST AREAS	SUGGESTED PATIENT USE
PINK	FTM065		1.8mm	26g	Low	Box of 100	Ward Blood Glucose Monitoring	Diabetes Clinics, Paediatrics, Geriatrics, Home Use
BLUE	FTM066		1.8mm	21g	Medium	Box of 100	Ward Blood Glucose Monitoring	Diabetes Clinics, most normal adult skins
BEIGE	FTM067		2.4mm	21g	Medium/High	Box of 100	Haemoglobin, cholesterol, blood typing, coagulation, & other diagnostic tests	Paediatrics and most normal adult skins
GREEN	FTM069		2.8mm	21g	High	Box of 100	Haemoglobin, cholesterol, blood typing, coagulation, & other diagnostic tests	thicker adult skins
STERIHEEL ORANGE: PREEMIE	FTM254		0.85 mm	1.75mm incision length	High	Box of 50	Special baby care and maternity units for heel puncture	Preemies
STERIHEEL YELLOW: NEWBORN	FTM253		1.0mm	2.5mm incision length	High	Box of 50	Special baby care and maternity units for heel puncture	Newborns

EQUIVALENTS TO VITALCARE LANCETS

AUTO TYPE	NPC CODE	HAEMOLANCE	SURGILANCE	UNISTIK	BD
PINK	FTM065	MICRO/LOW FLOW	YELLOW/GREY	COMFORT	GENIE – PURPLE
BLUE	FTM066	NORMAL FLOW	GREY	COMFORT/NORMAL	GENIE – PINK
BEIGE	FTM067	NORMAL/HIGH FLOW	ORANGE	NORMAL/EXTRA	GENIE – PINK
GREEN	FTM069	HIGH/MAX FLOW	PINK	EXTRA	GENIE – BLUE BLADE
STERIHEEL ORANGE	FTM254	PLUS PAEDIATRIC	TENDERFOOT PREEMIE		QUICKHEEL LAVENDER
STERIHEEL YELLOW	FTM253	PLUS PAEDIATRIC	TENDERFOOT NEWBORN		QUICKHEEL TEAL

Distributed by
Pat Coakley Medical Ltd
sales@patcoakleymedical.ie
www.patcoakleymedical.ie
phone: 01 2018070 / 2

