

TRACHEOSTOMY CARE

ADVADRAW (T) LAMINATED TRACHEOSTOMY DRESSING



LATEX FREE



Advadraw 'T' is a tri-laminate tracheostomy pad which has a low adherent perforated wound film contact layer, an absorbent centre action and a strike-through resistant layer to seal the dressing.

Advadraw 'T' is comprised of a new improved lighter fibrous material, with tests showing improved absorbency rates when compared to other manufacturers.

The slim contours of Advadraw 'T' avoids elevating the tracheostomy tube away from the neck, while the curved edges ensure maximum comfort to the patient.

Available in one convenient size to suit most intubation procedures, can be easily trimmed if needed.

- Cushioning effect on wound
- Low adherence -Reduced trauma on removal
- Conformable comfortable protection
- Suitable for low/medium levels of exudate
- Rapid uptake of exudates
- Available on FP10 Prescription
- Single patient use only
- UK Manufactured
- Latex-Free
- CE Marked



IMP\$CODE	DESCRIPTION	SIZE	DRUG\$ARIFF	NHS\$CODE	PACK\$QTY
CR/4416	Advadraw (T) Dressing	180 x 80mm	40994970	FDH615	20