

TRACHEOSTOMY TUBE CLEANING SWAB

The Insight Cleaning Swabs are perfect for cleaning all types of PVC and silicone tracheostomy or laryngectomy tubes.

Available in two head sizes, our Tube Cleaning Swabs are secure and suitable for popular sizes of tubes up to 6mm (small), 9mm (large) and 9mm (Extra Long Handle).

The swabs soft foam heads ensure damaging or scuffing the inside of the tube is avoided, whilst the white plastic handle is flexible for ease of use.



- Use with warm water to remove mucus residue in the tube
- Available in two sizes suitable for popular sizes of tube
- Available on FP10 Prescription
- Single patient use only
- Latex-Free
- CE Marked

IMP CODE	DESCRIPTION	SIZE	DRUG TARIFF	NHS CODE	PACK QTY
TCS/900	Tube Cleaning Swabs	Large	324-8689	FDH420	100 (10x10)
TCS/950	Tube Cleaning Swabs	Small	324-8671	FDH421	100 (10x10)
TCS/1000	Tube Cleaning Swabs	Extra Long	N/A	FDH611	100 (10x10)

IMP CODE	HEAD WIDTH	HEAD LENGTH	HANDLE LENGTH	MAXIMUM TUBE SIZE
TCS/900	13mm	26mm	155mm	9mm
TCS/950	8mm	22mm	155mm	6mm
TCS/1000	9mm	25mm	245mm	9mm