



All Silicone 2 way Foley Catheter

Description

All silicone indwelling Foley 2 way catheter, traditional balloon construction, hand finished to meet exact surgical and patient requirements. BS 1616.

Sizes

International	Male		Female		Paediatric
Colour Code					
F.G. Colour	5-10ml	30ml	5-10ml	30ml	3-5ml
8 Red					83 P
10 Yellow					83 P
12 White	84 M		86 F		
14 Green	84 M		86 F		
16 Orange	84 M	85 M	86 F	87 F	
18 Red	84 M	85 M	86 F	87 F	
20 Yellow	84 M	85 M	86 F	87 F	
22 Purple	84 M	85 M	86 F	87 F	
24 Blue	84 M	85 M	86 F	87 F	



Indication

For use when either short or medium term catheterisation is required.

Advantages

- Manufactured from high quality Silicone.
- Established design (Foley type).
- Silicone elastomer equates to low reactivity by body tissues, minimised urethral irritations and reduced catheterisation complications, such as trigonitis, urethritis and meatitis.
- Non-wetting surface facilitates drainage, is resistant to encrustation by urinary salts and is non-adherent, providing for greater patient comfort.
- Fewer catheter changes mean that infection caused by repeated removal is markedly reduced.
- Translucent with blue radio-opaque line.
- All our product range is latex and phthalate free.

Presentation

- Supplied sterile in individual double wrapped pouches, available in 1,3 and 10 per carton
- Perforated inner envelopes maximise protection against contamination.
- Also available as single unit.
- Published in NHS Drug Tariff.

Sterilization

- Sterilized for single use only. 3-5ml and 5-10ml balloon sterilized by gamma irradiation.
- 30ml balloon sterilized by ethylene oxide.