

Distributed by

Pat Coakley Medical Ltd

[sales@patcoakleymedical.ie](mailto:sales@patcoakleymedical.ie)

[www.patcoakleymedical.ie](http://www.patcoakleymedical.ie)

phone: 01 2018070 / 2



## Disinfection Tanks

Our range of EDPO polymeric tank-containers are designed to carry out the disinfection of medical equipment and pre-sterilisation treatment. These containers are made of high impact, chemical stable, non-toxic polystyrene and are compatible with the disinfectants and chemical agents used during these procedures with a maximum working temperature of 75°C.

The tank-containers consist of:

1. A body in the form of rectangular shape tank with an upper lip providing easy and convenient handling.
2. A perforated inner tray.
3. A press plate to allow full immersion into the solution for lighter equipment.
4. A lid to close the tank to prevent unwanted exposure to chemicals.

### Tank-Container Specifications:

<b>EDPO-1-01</b>	223mm x 149mm x 91mm	Utilised Volume 1 Litre
<b>EDPO-3-01</b>	315mm x 206mm x 125mm	Utilised Volume 3 Litre
<b>EDPO-5-01</b>	394mm x 260mm x 156mm	Utilised Volume 5 Litre
<b>EDPO-10-01</b>	496mm x 328mm x 195mm	Utilised Volume 10 Litre
<b>EDPO-10D-01</b>	870mm x 235mm x 160mm	Utilised Volume 10 Litre

