



## Laserguard Ocular Range

### Description

A range of laser resistant ocular contact shields for use in periorbital dermatological and cosmetic laser treatments.

### LGE 1 Y

A yellow eye shield that has been tested and found safe for use with PULSE DYE LASER (YELLOW/ORANGE) LIGHT at 585 to 600nm.

### LGE 1 R

A red eye shield that has been tested and found safe for use with RUBY LASER (RED) LIGHT at 694.7nm and with ALEXANDRITE LASER LIGHT at 752nm.

### LGE 1 G

A green eye shield that has been tested and found safe for use with Nd:YAG LASER (GREEN) LIGHT at 532nm.

### Presentation

- Supplied sterile in individual double wrapped pouches, 25 pouches per carton.

### Sterilization

- Sterilized by gamma irradiation, single use only.

### LGE 2

- A stainless steel ocular contact shield that has been tested and found safe for use with ERBIUM YAG AND CO<sub>2</sub> ULTRA PULSE LASER (INFRA RED) LIGHT at 1.064 microns to 10.64 microns.

### Presentation

- Supplied non-sterile in single wrapped pouch, two shields per pouch.
- The ocular contact shield can be cleaned and re sterilized using recognised validated hospital instrument cycles.

### LGE 2A

- Silicone applicator designed for use with LGE2 stainless steel, ocular contact shield.

### Presentation

- Supplied non-sterile in individual single pouch. Single item.

### Sterilization

- RENEW IF DIRECT STRIKE RECEIVED.  
NOTE: The fitting of any eye shield requires an anaesthetic, either local or general.
- For sterilization instructions see page 41

### Latex & Phthalate Free

- All our product range is latex and phthalate free.





---

## Laser Resistant Bite-Guard

---

### Description

A blue bite guard that has been tested and found safe for use with Erbium Yag laser (Infra Red) light at 2.94 microns.

### Sizes

Code	Size
LBG1	One size only

### Indication

A Laser Resistant Bite-Guard to protect dental enamel during Erbium Yag laser procedures to the face.

### Advantages

- Found safe for use with Erbium and Yag laser (Infra Red) light at 2.94.
- Manufactured from high quality medical grade silicone.
- Radio opaque.
- All our product range is latex and phthalate free.

### Presentation

- Supplied sterile in individual double wrapped pouches, 20 pouches per carton.

### Sterilization

- Sterilized by gamma irradiation, Single use only.

