



M. C. Johnson Company, Inc.



Distributed by

Pat Coakley Medical Ltd

sales@patcoakleymedical.ie

www.patcoakleymedical.ie

phone: 01 2018070 / 2



WOUND • SECURE[®]

ADHESIVE FOAM DRESSINGS

ALGINATE DRESSINGS

BORDERED GAUZE DRESSINGS

HYDROCOLLOID DRESSINGS w/HYDROWEB 18™

NON-ADHESIVE FOAM DRESSINGS

RETENTION TAPE

TRANSPARENT DRESSINGS

WATERPROOF DRESSINGS

EAR LOOP MASK w/SHIELD

FULL-FACE SHIELD

CATH • SECURE[®]


CATH-SECURE[®]-2

CATH-SECURE[®]-3

CATH-SECURE DUAL TAB[®]-4

CATH-SECURE PLUS[®]-6

NG SECURE[®]



The NG SECURE® is made using a hypoallergenic tape and features a moisture-resistant foam cushion that creates a barrier between the skin and the tubing. By creating this barrier, the foam cushion on NG SECURE® keeps pressure off of the patient's nostrils and lips which may help to prevent pressure sores and discomfort. Taping to the oily nose area, conventional tube holders fall off within hours. The device has a resealable tape flap for securing the tube onto the top of the foam cushion positioned under the nose. This innovative design permits the tubing to be readjusted or repositioned as necessary without retaping and allows for the usage of simultaneous multiple tubing's.

Key Features:

- Conforms to OSHA Standard 1910.1030
- Latex Free
- Comfortable
- Easy to Apply
- Readjustable
- Versatile
- Securely holds feeding and/or aspiration tubes in place
- Gentle on delicate skin
- Resists humidified air of oxygen mask
- Adhesive is not affected by oily skin
- Leaves no sticky residue on tube or nose after removal

Compares to:

- Stabilock - Dale
- NG Tube Holder - Mormac TubeGuard
- Nasogastric Tube Holder, Osfa – Dale



LATEX FREE • DEHP-FREE • HYPOALLERGENIC

Size

5" x 1" base, 1.875" x .375" holder

Package

50 BOX / 8 CASE


Product

9304-1

MANUFACTURED FOR:
M. C. JOHNSON CO., INC.
8801 Business Park Drive
Suite 201
Ft. Myers, FL. 33912

NG•SECURE®
M.C. JOHNSON®
M. C. Johnson Company, Inc.
www.mcjohnson.com

PH: (800) 553-8483
PH: (239) 433-0619
FAX: (239) 433-1498
FAX: (800) 643-0343




The TUBE-SECURE *for Kids*™ is a multi-purpose tube holder. It provides an efficient, standardized system for securing pediatric medical tubing. The device is hypoallergenic, DEHP-FREE and Latex Free. The TUBE-SECURE *for Kids*™ lowers healthcare costs and increases nursing productivity while providing quality care for patients and meeting the healthcare industry standard of care. The TUBE-SECURE *for Kids*™ can be purchased through your preferred medical distributor.

Key Features:

- Improves patient comfort and mobility
- Eliminates the need for constant re-taping
- Permits repeated access to tubing without restricting the flow
- Quick and easy to apply
- Allows patient movement without fear of pulling catheter
- Flexible and lightweight

Compares to:

- Bard – Pediatric Statlock®



LATEX FREE • DEHP-FREE • HYPOALLERGENIC

Securement Areas (not limited to):

- Indwelling urinary catheters
- Gastrostomy tubes (G-Tubes)
- Jejunostomy tubes

Size	Package	HCPSCS	Product #
1.87 X .87" Base, 1.31" x .25" Tab	50 BOX / 8 CASE	A4333 or A5200	7334-1


MANUFACTURED FOR:
M. C. JOHNSON CO., INC.
8801 Business Park Drive
Suite 201
Ft. Myers, FL. 33912

TUBE SECURE *for Kids*™



M. C. Johnson Company, Inc.
www.mcjohnson.com

PH: (800) 553-8483
PH: (239) 433-0619
FAX: (239) 433-1498
FAX: (800) 643-0343




The CATH•SECURE Plus™ is made using a breathable, water-resistant, hypoallergenic material that features a butterfly design. The device above is DEHP-Free and Latex Free. This product line improves overall clinical outcomes, lowers healthcare costs, increases nursing productivity all while providing a better quality of care for patients and meeting the healthcare industry standard of care.

Key Features:

- Improves patient comfort and mobility
- Eliminates the need for constant re-taping
- Permits repeated access to tubing without restricting the flow
- Quick and easy to apply
- Allows patient movement without fear of pulling catheter
- Flexible and lightweight
- JCAHO, OSHA, and CDC compliant

Compares to:

- Venetec – Statlock®



LATEX FREE • DEHP-FREE • HYPOALLERGENIC

Securement Areas (not limited to):

- Indwelling urinary catheters
- Gastrostomy tubes (G-Tubes)
- Jejunostomy tubes

Size	Package	HCPCS	Product #
2.5" Long Tab & Butterfly Base	50 BOX / 8 CASE	A4333 / A5200	5445-6

MANUFACTURED FOR:
M. C. JOHNSON CO., INC.
8801 Business Park Drive
Suite 201
Ft. Myers, FL. 33912


CATH•SECURE Plus™



M. C. Johnson Company, Inc.
www.mcjohnson.com

PH: (800) 553-8483
PH: (239) 433-0619
FAX: (239) 433-1498
FAX: (800) 643-0343

CATH•SECURE DUAL TAB®




The CATH-SECURE DUAL TAB® is designed to extend from the center, wrap around the tubing, and anchor it tightly to the tape using a press-on hook-and-loop fastening system. The CATH-SECURE DUAL TAB® improves patient comfort and mobility, is lightweight and flexible, and most of all is JCAHO, OSHA and CDC compliant.

Key Features:

- Eliminates the need for constant re-taping
- Permits repeated access to tubing without restricting the flow
- Quick and easy to apply
- Allows patient movement without fear of pulling catheter
- Quick and easy to remove
- Cost effective

Wrap around tubing and press onto hook and loop fastener



Wrap each dual tab strip around and under tubing

LATEX FREE • DEHP-FREE • HYPOALLERGENIC

Securement Areas (not limited to):

- Peritoneal Dialysis exit site tubing's
- Leg bag tubing
- Condom catheters
- Prostatectomy tubes

Size	Package	HCPCS	Product #
(2-Dual) 3" Wrap Tabs x 3" Wide Base	50 BOX / 8 CASE	A4333 / A5200	5445-4


MANUFACTURED FOR:
M. C. JOHNSON CO., INC.
8801 Business Park Drive
Suite 201
Ft. Myers, FL. 33912

CATH•SECURE DUAL TAB®

M.C. JOHNSON
M. C. Johnson Company, Inc.
www.mcjohnson.com

PH: (800) 553-8483
PH: (239) 433-0619
FAX: (239) 433-1498
FAX: (800) 643-0343

Rev. 2




The CATH-SECURE® For Kids is a multi-purpose tube holder. It provides an efficient, standardized system for securing pediatric catheters and medical tubing. The device is hypoallergenic, DEHP-FREE and Latex Free. The CATH-SECURE® For Kids lowers healthcare costs and increases nursing productivity while providing quality care for patients and meeting the healthcare industry standard of care.

Key Features:

- Improves patient comfort and mobility
- Eliminates the need for constant re-taping
- Permits repeated access to tubing without restricting the flow
- Quick and easy to apply
- Allows patient movement without fear of pulling catheter
- Flexible and lightweight

Compares to:

- Bard – Pediatric Statlock®



LATEX FREE • DEHP-FREE • HYPOALLERGENIC

Securement Areas (not limited to):

- Indwelling urinary catheters
- Gastrostomy tubes (G-Tubes)
- Jejunostomy tubes

Size	Package	HCPSC	Product #
1¼" Long Tab & 1½" x ¾" Wide Base	50 BOX / 8 CASE	A4333 or A5200	7334-2


MANUFACTURED FOR:
M. C. JOHNSON CO., INC.
8801 Business Park Drive
Suite 201
Ft. Myers, FL. 33912

CATH•SECURE® for Kids



M. C. Johnson Company, Inc.
www.mcjohnson.com

PH: (800) 553-8483
PH: (239) 433-0619
FAX: (239) 433-1498
FAX: (800) 643-0343




CATH-SECURE™ is the original multi-purpose tube holder. It provides an efficient, standardized system for securing catheters and medical tubing. The device is easy to adjust and may remain in place for approximately one week without altering skin integrity. CATH-SECURE™ is made with alcohol-soluble adhesive.

Key Features:

- Holds different sizes of tubing, as well as multiple tubes
- Improves patient's quality of life
- Quick and easy to remove
- Cost effective
- Increases nursing productivity
- OSHA, and CDC compliant

Compares to:

- Venetec – Statlock®



LATEX FREE • DEHP-FREE • HYPOALLERGENIC CE 2

Securement Areas (not limited to):

- Subclavian line ports
- EKG lead wires
- Dialysis tubing

Size	Package	HCPCS	Product #
2.5" Long Tab & 3" x 1.75" Wide Base	50 BOX / 8 CASE	A4333 / A5200	5445-2
2.75" Long Tab & 3" x 1.25" Wide Base	50 BOX / 8 CASE	A4333 / A5200	5445-3

MANUFACTURED FOR:
 M. C. JOHNSON CO., INC.
 8801 Business Park Drive
 Suite 201
 Ft. Myers, FL. 33912


CATH•SECURE™

 M. C. Johnson Company, Inc.
www.mcjohnson.com


PH: (800) 553-8483
 PH: (239) 433-0619
 FAX: (239) 433-1498
 FAX: (800) 643-0343

CATH•SECURE™

The Original Multi-Purpose Tube Holder



Don't Settle for an Imitation



Features

- CATH•SECURE™ provides an efficient standardized system for securing catheters and other medical tubing.
- CATH•SECURE™ is a breathable, hypoallergenic material that features a butterfly design.
- CATH•SECURE™ may comfortably remain in place for several days without effecting skin integrity.
- CATH•SECURE™ is made with an alcohol-soluble adhesive.
- CATH•SECURE™ is easy to adjust and remains securely in place.
- CATH•SECURE™ holds different sizes of tubing, as well as multiple tubes.

Benefits

- CATH•SECURE™ may reduce the rate of urinary tract infections by up to 45%.
- CATH•SECURE™ improves patient comfort and mobility.
- CATH•SECURE™ increases nursing productivity.
- CATH•SECURE™ improves the patient's quality of life.
- CATH•SECURE™ eliminates the need for constant re-taping.
- CATH•SECURE™ permits repeated access to tubing without restricting the flow.
- CATH•SECURE™ is JCAHO, OSHA and CDC Compliant.

Product #	Description	Package	HCPCS
5445-2	CATH-SECURE™ with 2.5" Tab	50/Box	A4333
5445-3	CATH-SECURE™ with 2.75" Tab	50/Box	A4333
5445-4	CATH-SECURE DUAL TAB®	50/Box	A4333
5445-5	Lite & Breathable CATH-SECURE™	50/Box	A4333
5445-6	CATH-SECURE Plus™	50/Box	A4333



LATEX FREE • DEHP-FREE • HYPOALLERGENIC • F315 COMPLIANT

MANUFACTURED FOR:
M. C. JOHNSON CO., INC.
8801 Business Park Drive
Suite 201
Ft. Myers, FL. 33912

CATH•SECURE™

M. C. JOHNSON®

M. C. Johnson Company, Inc.

www.mcjohnson.com

PH: (800) 553-8483

PH: (239) 433-0619

FAX: (239) 433-1498

FAX: (800) 643-0343